



Be Right™



USB + Power Module, 2100Q and DR 1900

Proizvod #: LZV813.99.00001

HRK Cijena: Kontaktirajte nas

Otprema unutar 1 tjedna

The USB + Power Module is an optional accessory for the 2100Q & 2100Q IS portable turbidimeters, and the portable spectrophotometer DR1900. This module adds further functionality to the instrument allowing the user to:

Transfer data to a PC in XML format

Transfer data to the Citizen PD-24 printer (2960100)

Operate the instrument on line power (100-240Vac, 50/60 HZ) and Charge NiMH batteries (2971304)

Perform firmware updates

Quickly enter Sample and User IDs with the optional barcode scanner (LZV566)

Specifikacije

Output Voltage:	9V
Uključuje :	Includes module, module cover, power supply with four plug adapters, standard USB cable with Mini-USB connector, special USB cable adapter, and manual.
Ulazni napon:	100 - 240 VAC, 50/60 Hz