




pH electrode, diameter 12 mm, 10 m fixed moulded cable, pHRET

Proizvod #:

LZX476

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Electrode for the determination of the pH value in waste water and pressurised pipes. Optionally with integrated temperature sensor. Suitable for use with the 12 mm inline armature.

Specifikacije

Dijafragma:	Porous PTFE, double ring diaphragm
Dubina uranjanja:	205 mm
Duljina:	225 mm
Duljina kabla:	10 m fixed cable
Elektrolit:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garancija:	6 mjeseci
Materijal:	Stainless steel
Mjerni rang:	0 - 14 pH
Radna temperatura:	-5 - 135 °C
Rang tlaka:	15 bar
Reference system:	Ag/AgCl
Temperaturni senzor:	None