



pHRETPT100 pH sensor, 12 mm, retractable

Proizvod #:	LZX477
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

With integrated temperature sensor PT100 and 10 m cable. Suitable for pressure lines.

Specifikacije

Dijafragma:	Porous PTFE, double ring diaphragm
Dubina uranjanja:	205 mm
Duljina:	225 mm
Duljina kabla:	10 m fixed cable
Elektrolit:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garancija:	6 mjeseci
Materijal:	Stainless steel
Mjerni rang:	0 - 14 pH
Radna temperatura:	-5 - 135 °C
Rang tlaka:	15 bar
Reference system:	Ag/AgCl
Temperaturni senzor:	PT100