



LTLCONT pH Electrode with Top68 Connector

Proizvod #: LZX536
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Low temperature glass and electrolyte, single reference system, electrode for ion-weak and/or cold water (10-100 μ S/cm), without temperature sensor

Specifikacije

Dijafragma:	Ring diaphragm, porous PTFE
Dubina uranjanja:	105 mm
Duljina:	140 mm
Elektrolit:	KCl, saturated KCl/AgCl (crystal)
Garancija:	6 mjeseci
Mjerni rang:	2 - 11 pH
Navoj senzora:	Pg 13,5
Radna temperatura:	-20 - 50 °C
Rang tlaka:	3,5 bar
Reference system:	Ag/AgCl
Spajanje kablova:	TOP68 connector
Temperaturni senzor:	None