



SOTA10 pH electrode, 10 m fixed cable, without T-sensor

Proizvod #: LZX587
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Standard glass and reference system

Specifikacije

Dijafragma:	Ring diaphragm, porous PTFE, double
Dubina uranjanja:	105 mm
Duljina:	140 mm
Duljina kabla:	10 m fixed cable
Elektrolit:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garancija:	6 mjeseci
Mjerni rang:	0 - 14 pH
Radna temperatura:	0 - 80 °C
Rang tlaka:	3,5 bar