




pH glass electrode, refillable, pure water applications

Proizvod #:	LZY025
EUR Cijena:	Kontaktirajte nas
Dostupno	

pH electrode for use in ultra-pure water

Specially developed for use in ultra-pure water

Refillable glass electrodes for special requirements and high-precision measurements

Improved durability through triple ceramic diaphragm

Easy to install in stainless steel or polypropylene flow fittings

Specifikacije

Dijafragma:	Triple ceramic
Duljina:	120 mm
Elektrolit:	KCl liquid
Garancija:	6 mjeseci
Masa:	0,15 kg
Materijal osovine elektrode:	Glass
Mjerni rang:	0 - 12 pH
Radna temperatura:	-30 - 80 °C
Reference system:	Ag/AgCl
Spajanje kablova:	Plug head S7
Temperaturni senzor:	None
Vrijeme odziva T98%:	< 30 s