



pH Electrode for 8362 Ultrapure Water Sensor

Proizvod #: 08362=A=2000
EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

Part of analogue pH sensor 8362, specially designed for pure water pH measurement.

Specifikacije

Duljina:	120 mm
Garancija:	6 mjeseci
Konektori:	S7 Connector
Materijal:	Sensor: Glass
Mjerni rang:	0 - 14 pH
Navoj senzora:	Pg 13,5
Posebne napomene:	pure water pH electrode
Primjena:	Ultrapure water
Referentna elektroda:	None
Spajanje kabla:	S7 connector
Temperaturni senzor:	PT100