



## Recharge kit of 4 pre-filled cartridges with premounted 2935A membranes for A1100, GA2400 & GA2800 EX oxygen sensors

**Proizvod #:**

2935A-A

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna

Opasno po zdravlje



Includes o-rings, cotton filters and cleaning tools.

---

### Specifikacije

Kompatibilni instrumenti:

A1100, GA2400 & GA2800 EX oxygen sensors

Membrana:

2935A membranes