



Standard za nitrat, 200 mg/L NO₃ (45,2 mg/L NO₃-N), 500 mL

Proizvod #:

LCW827

EUR Cijena:
Dostupno

Kontaktirajte nas

Standardna otopina za online određivanje nitrata

Specifikacije

Koncentracija:	200 mg/L NO ₃ (45,2 mg/L NO ₃ -N)
Parametar:	Nitrat
Platforma:	Otopina
Uređaj:	Nitratax
Volumen/ veličina paketa:	500 mL