



Beaker, polypropylene, low form, 50 mL

Proizvod #:	108041
EUR Cijena:	Kontaktirajte nas
Dostupno	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Jedinica:	1 each
Kapacitet:	50 mL
Materijal:	Polypropylene, Low Form
Materijal tijela mjerača:	Polypropylene
Promjer, vanjski:	51.1 mm with spout
Volumen/ veličina paketa:	1 each