



## Nitrate Standard Solution, 1000 mg/L, 500 mL

**Proizvod #:**

**1279249**

**Opasno po zdravlje**

EUR Cijena:

Kontaktirajte nas

Dostupno



### Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Nitrate, 1000 mg/L  $\pm 10$  mg/L  $\text{NO}_3\text{-N}$  (NIST). 500 mL bottle

#### Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifikacije

Čuvanje uređaja:	10 - 25 °C (protect from light)
Koncentracija:	4430 $\pm$ 44 mg/L $\text{NO}_3$
Opis:	Nitrogen-Nitrate Standard Solution as $\text{NO}_3\text{-N}$ (NIST)
Package settings:	Bottle
Parametar:	Nitrate
Platforma:	Otopina
Posebne napomene:	Observed shelf life after opening, following recommended storage and use: 3 months
Volumen/ veličina paketa:	500 mL