



Be Right™



Copper Standard Solution, 10 mg/L Cu (NIST), 100 mL MDB, for APA6000

Proizvod #: 12932
EUR Cijena: Kontaktirajte nas
Dostupno

Specifikacije

Čuvanje uređaja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Koncentracija:	10 mg/L
Opis:	Copper Standard Solution as Cu (NIST)
Parametar:	Hardness
Uređaj:	APA6000
Volumen:	100 mL MDB
Volumen/ veličina paketa:	100 mL