



Volatile Acid Standard Solution, 62,500 mg/L as Acetic Acid, pk/16 - 10 mL Voluette Ampules

Proizvod #:	1427010
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Acetate Standard Solution. For accuracy checks in volatile acids determinations.

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specifikacije

Koncentracija:	62500 mg/L
Platforma:	Otopina
Volumen/ veličina paketa:	10 mL/ pack of 16