



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 250 mL, pk/12

Proizvod #: 1472456
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikacije

Body Material:	Plastic
Kapacitet:	250 mL
Materijal:	PP
Package size UOM:	12pk
Volumen/ veličina paketa:	12 /pk