



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

Proizvod #:	1472459
EUR Cijena:	Kontaktirajte nas
Otprema unutar 2 tjedna	

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikacije

Body Material:	Plastic
Kapacitet:	500 mL
Materijal:	PP
Package size UOM:	12pk
Volumen/ veličina paketa:	12 /pk