



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 1 L, 6/pk

Proizvod #: 1472463
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikacije

Body Material:	Plastic
Kapacitet:	1000 mL
Materijal:	PP
Package size UOM:	6pk
Volumen/ veličina paketa:	6 /pk