



Cyanuric acid test kit, model CY-3

Proizvod #: 185102
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Broj testova:	100
Mjerni rang:	20 - 100 mg/L
Model:	CY-3
Naziv metode:	Turbidimetry
Parametar:	Cyanuric Acid
Primjena:	Pool and Spa
Smallest increment:	20
Test Kit Model:	CY-3
Vrsta kućišta:	D