



Cylinder, graduated, mixing, glass, 25 mL ± 0.3 mL, 0.5 mL divisions, glass stopper #13

Proizvod #: 189640
EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specifikacije

Certificiranja:	Class B
Kapacitet:	25 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Contain
Volumen/ veličina paketa:	1 each
Vrsta:	Jednometrički, miješanje