



EDTA standard solution, 0.0200N, 4L

Proizvod #: 20556
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

For hardness determination by buret titration. Hach methods 8222 (calcium) and 8226 (total). 0.0200 N (0.0100 M) (N/50) 4 liter bottle

Specifikacije

Mjerni rang:	0 - 5,000 mg/L CaCO_3
Naziv metode:	EDTA Titration
Platforma:	Otopina
Volumen/ veličina paketa:	4

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)