



APA6000 Hardness Standard Solution, 0.50 mg/L CaCO₃ (NIST), 946 mL

Proizvod #:	2058016
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Specifikacije

Čuvanje uređaja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Koncentracija:	0,50 mg/L
Opis:	Hardness Standard Solution (NIST)
Parametar:	Hardness
Uređaj:	APA6000
Volumen:	946 mL