



Bottle, BOD, black, 300 mL

Proizvod #:	2086047
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

Specifikacije

Body Material:	Glass
Kapacitet:	300 mL
Materijal:	Borosilicate Glass
Package size UOM:	1pk
Volumen/ veličina paketa:	1 each