



Bottle, Sampling, Low Density, Polyethylene, 250 mL, 12/pk

Proizvod #:	2087076
EUR Cijena:	Kontaktirajte nas
Otprema unutar 2 tjedna	

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Specifikacije

Body Material:	Low Density Plastic
Kapacitet:	250 mL
Materijal:	LDPE
Package size UOM:	12pk
Volumen/ veličina paketa:	12 /pk