



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 1000 mL

Proizvod #:	2089653
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specifikacije

Jedinica:	1 each
Kapacitet:	1000 mL
Materijal:	Glass Screw Cap