



Flask, Erlenmeyer, polymethylpentene capacity 125 mL

Proizvod #:	2089843
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specifikacije

Jedinica:	1 each
Kapacitet:	125 mL
Materijal:	Polymethylpentene
Vrsta:	Erlenmeyer