



Be Right™



Flask, Erlenmeyer, polymethylpentene capacity 250 mL

Proizvod #:	2089846
EUR Cijena:	Kontaktirajte nas
Dostupno	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specifikacije

Jedinica:	1 each
Kapacitet:	250 mL
Materijal:	Polymethylpentene
Vrsta:	Erlenmeyer