



Pipeta, serološka, sterilna, od polistirena, 5 mL, plava, komplet od 200 komada

Proizvod #:	2092667
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Sterilna plastična pipeta za jednokratnu upotrebu. Pojedinačno pakirano. Oznake u boji. Kalibrirano na izljev. U skladu s ASTM E934.

Graduations are calibrated for accurate dispensing to within $\pm 2\%$

Pipets meet BST/TSE requirements; latex-free

DNase/RNase-free, non-pyrogenic, non-cytotoxic, non-haemolytic

Specifikacije

Graduations:	0.10 mL
Kapacitet:	5 mL
Negative Graduations:	2.0 mL
Volumen/ veličina paketa:	200 u kompletu