



Reagent powder pillows, CuVer 1 copper, 10 mL, pk/1000

Proizvod #:

2105828

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For copper determination by the bichinchoninate method. Range: up to 5.00 mg/L Cu. EPA compliant for wastewater reporting - Hach Method 8506. Sample Size: 10 mL. Pack of 1000 powder pillows.

Specifikacije

Broj metode:	8506
Broj testova:	1000
EPA compliant:	Da
Mjerni rang:	0,04 - 5,00 mg/L Cu
Naziv metode:	Bichinchoninate
Parametar:	Copper
Platforma:	Praškasti reagensi PP
Uređaj:	All except DR/3000, DR/2000, DR/850, DR/820
Volumen/ veličina paketa:	1000 /pk