



## Standard solution, manganese, 25 mg/L as Mn (NIST), pk/20 - 2 mL PourRite ampules

**Proizvod #:**

**2112820**

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



For accuracy checks in manganese determinations

### Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

---

## Specifikacije

Koncentracija:	25 mg/L
Platforma:	Otopina
Volument/ veličina paketa:	mL/20