



Be Right™



Centrifuge Tube, Conical Bottom, 10/pk

Proizvod #:	2115839
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

Specifikacije

Duljina:	120 mm
Kapacitet:	15 mL
Volumen/ veličina paketa:	10 /pk