



Nitrogen-Ammonia Standard Solution, 150 mg/L NH₃-N, pk/16 - 10 mL Voluette Ampules

Proizvod #: 2128410
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

For accuracy checks in TKN determinations

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specifikacije

Koncentracija:	150 mg/L
Opis:	Ammonia Standard Solution as NH ₃ -N
Platforma:	Otopina
Volumen/ veličina paketa:	mL/16