



Bag, sterile, WhirlPak, 207 mL, 500/pk

Proizvod #: 2233100
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Approved by U.S. Environmental Protection Agency (USEPA). Disposable, transparent polyethylene bags. Gas-sterilized, pull-tab openings. Can be marked with waterproof felt pen or marker. Sure "whirl" closure. Nasco.

Specifikacije

Kapacitet:	207 mL
Volumen/ veličina paketa:	500