



Tube, culture, 16 x 100 mm, pk/6

Proizvod #:	2275806
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

CAPS SOLD SEPARATELY product number 2241106 (6 pieces).

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Specifikacije

Duljina:	100 mm
Kapacitet:	Capacity to Neck 11 mL
Promjer:	16 mm
Volumen/ veličina paketa:	6 pk