



Agar R2A, dehydrated, 500 g

Proizvod #:	2281034
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Low nutrient medium for enumeration and cultivation of bacteria from potable water. Used for pour or spread plate and MF methods; described in the 16th edition of Standard Methods. Recommended for predicting nitrification more accurately in potable systems using chloramines as a secondary disinfectant.

Specifikacije

Media Type:	R2A Agar
Volumen/ veličina paketa:	500 g