



Conductivity Standard Solution, 180 $\mu\text{S}/\text{cm}$, NaCl, 1000 mL

Proizvod #: 2307553
EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 180 $\mu\text{S}/\text{cm} \pm 10 \mu\text{S}/\text{cm}$ at 25°C, 85.47 mg/L as Sodium Chloride (NaCl), (1 $\mu\text{S}/\text{cm} \sim 0.5 \text{ mg/L TDS}$), 1000 mL bottle.

Specifikacije

Čuvanje uređaja:	10 - 25 °C (protect from light)
Formula:	NaCl
Koncentracija:	180 $\mu\text{S}/\text{cm}$ at 25°C $\pm 10 \mu\text{S}/\text{cm}$ at 25°C
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parametar:	Conductivity
Platforma:	Otopina
Posebne napomene:	Observed shelf life after opening, following recommended storage and use: 3 months
Volument/ veličina paketa:	1000 mL