



Butterfield's buffered phosphate, 450 mL, pk/10

Proizvod #:	2319110
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Butterfields Buffered Phosphate, 450 mL bottle. pk/10.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Specifikacije

Jedinica:	10/pk
Volumen:	450 mL
Volumen/ veličina paketa:	10/pk