



Be Right™



Test kit, oxygen scavenger, color disc, model OS-1, 100 tests

Proizvod #:

2349500

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 100 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Broj testova:	50
Mjerni rang:	0 - 52 µg/L DEHA (and other Oxygen Scavengers)
Mjerni rang (2):	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Model:	OS-1
Naziv metode:	Color disc/Iron reduction
Parametar:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platforma:	Disk sa bojama-Color Disc
Što je u pretincu:	173200 Color comparator box 1419705 Glass dropper w/ 0.5 & 1.0 mL marks, pk/5 2128800 Viewing tube insert 2412200 Adapter, long path viewing 4660004 Viewing tubes with caps, plastic, 4 pieces 1704200 Boca, čaša za miješanje/raspodjelu, 25 mL 9265900 Color Disc Oxygen Scavenger, 0 - 80 units

2167969 Jastučíci s DEHA reagensom 1 u prahu, 3 - 450 µg/L DEHA

2168042 Otopina DEHA reagens 2 (100 mL)

4660400 MultiTest kit case (142 x 312 x 185 mm), blue polypropylene

Test Kit Model:

OS-1

Vrsta kućišta:

R

Što je u pretincu

173200 Color comparator box
1419705 Glass dropper w/ 0.5 & 1.0 mL marks, pk/52
128800 Viewing tube insert
2412200 Adapter, long path viewing
4660004 Viewing tubes with caps, plastic, 4 pieces
1704200 Boca, čaša za miješanje/raspodjelu, 25 mL
9265900 Color Disc Oxygen Scavenger, 0 - 80 units
2167969 Jastučíci s DEHA reagensom 1 u prahu, 3 - 450 µg/L DEHA
2168042 Otopina DEHA reagens 2 (100 mL)
4660400 MultiTest kit case (142 x 312 x 185 mm), blue polypropylene