



Filter, acrodisc, pore size 0.2 µm

Proizvod #:	2415600
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Specifikacije

Filter Type:	Filter, acrodisc
Materijal:	Polysulfone
Packing:	NA
Promjer:	NA
Veličina pora:	0.2 µm