



Reagent set, copper, free & total

Proizvod #:

2439200

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



For determination of free and complexed (chelated) Copper by the Binchoninate method. Range: up to 5.00 mg/L. Set includes: Free Copper Reagent Powder Pillows and Hydrosulfite Powder Pillows for 100 tests.

Specifikacije

Broj testova:	100
Mjerni rang:	0,04 - 5,00 mg/L Cu
Naziv metode:	Bicinchoninate
Platforma:	Praškasti reagensi PP
Uređaj:	DR 6000, DR 3900, DR 1900, DR 900, PCII