



Reagent set, chloride, LR, digital titrator

Proizvod #:

2444100

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



For determination of Chloride by argentometric Digital Titrator titration. Range: 10 - 100 mg/L Cl⁻. Reagent Set contains Chloride 2 Indicator Powder Pillows and Silver Nitrate Digital Titrator Cartridges (2) 0.2256 N. Approximately 100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Broj metode:	8207
Broj testova:	100
Mjerni rang:	10 - 100 mg/L Cl ⁻
Naziv metode:	Silver Nitrate Titration