



Hardness (Calcium & Total Sequential) LR Reagent Set, Digital Titrator

Proizvod #:

2448600

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



Optimize Operational Efficiency with Easy-to-Use Water Hardness Test Kits

Water hardness test kits are essential tools in industrial and commercial environments for monitoring water quality and protecting equipment. These kits provide a convenient way to measure water hardness both in the lab and in the field, ensuring that hardness levels remain within acceptable ranges. By preventing scaling in pipes and boilers, they help optimize the performance of soaps and detergents, contributing to overall operational efficiency.

Our water hardness test kits are trusted by professionals in water treatment, food and beverage, textile, manufacturing, and other industries.

For sequential determination of Calcium and Total Hardness by EDTA Digital Titrator titration. Hach method 8329. Range: 10-160 mg/L CaCO₃. Set contains Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), Sulfuric Acid (5.25 N) and 0.0800 M EDTA Digital Titrator Cartridge. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Broj metode:	8329
Broj testova:	40 - 100
Mjerni rang:	10 - 160 mg/L CaCO ₃
Naziv metode:	EDTA Titration
Parametar:	Hardness (Calcium and Total Sequential)
Platforma:	Digitalni Titrator
Posebna značajka:	Reagent set, LR, Digital Titrator

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)