



Test tube, rack, polyethylene, 25 mm (O.D.), 24 holes

Proizvod #:	2497903
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifikacije

Aranžman testova:	3 x 8
Promjer:	25 mm
Rupe za testove:	24
Volumen/ veličina paketa:	1 each