



Test Tube, rack, polyethylene, 20 mm (O.D.), 20 holes

Proizvod #:	2497912
EUR Cijena:	Kontaktirajte nas
Dostupno	

Twenty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may Vary.

Specifikacije

Aranžman testova:	4 x 5
Promjer:	20 mm
Rupe za testove:	20
Volumen/ veličina paketa:	1 each