



## Sample container, 120 mL, pk/250

|                         |                   |
|-------------------------|-------------------|
| Proizvod #:             | 2552212           |
| EUR Cijena:             | Kontaktirajte nas |
| Otprema unutar 2 tjedna |                   |

Plastic, leak-proof. Corning Snap Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

---

### Specifikacije

|                           |                 |
|---------------------------|-----------------|
| Kapacitet:                | 120 mL - 250/pk |
| Promjer, vanjski:         | 45 mm (approx)  |
| Volumen/ veličina paketa: | 250 /pk         |