



Sample vial, 120 mL w/out thiosulfate, 200/pk

Proizvod #:	2634002
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Specifikacije

Duljina:	93 mm
Kapacitet:	Fill line at 100 mL ($\pm$ 2.5 mL)
Posebna značajka:	Sample Vial, without Thiosulfate
Promjer:	43 mm
Volumen/ veličina paketa:	200/pk