



Thermometer, Enviro-Safe, Non-Mercury, -10 - 260 °C

Proizvod #:	2635701
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Specifikacije

Duljina:	405 mm
Podjela termometra:	1 °C
Temperaturni rang:	-10 °C - 260 °C
Tip termometra:	Non-Mercury
Volumen/ veličina paketa:	1 each