



Cylinder, Graduated, 25 mL, Certified

Proizvod #: 2636240
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Specifikacije

Certificiranja:	Class A
Kapacitet:	25 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Deliver
Visina:	192 mm
Volumen/ veličina paketa:	1 each
Vrsta:	Jednometrički