



## Cylinder, Graduated, 500 mL, Certified

**Proizvod #:** 2636249  
**EUR Cijena:** Kontaktirajte nas  
**Nije prikazan datum isporuke**

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes  $\geq 25$  mL.

---

### Specifikacije

|                           |               |
|---------------------------|---------------|
| Certificiranja:           | Class A       |
| Kapacitet:                | 500 mL        |
| Materijal tijela mjerača: | Glass         |
| Mjerna metoda:            | To Deliver    |
| Visina:                   | 390 mm        |
| Volumen/ veličina paketa: | 1 each        |
| Vrsta:                    | Jednometrički |