



Cylinder, graduated, mixing, certified, 100 mL

Proizvod #:	2636342
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Specifikacije

Certificiranja:	Class A
Kapacitet:	100 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Contain
Volumen/ veličina paketa:	1 each
Vrsta:	Jednometrički, miješanje