



Clamp, Extension Thermometer

Proizvod #:	2656200
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Replacement clamp for use with the Cold Vapor Mercury Test. Hold thermometers and other glass apparatus firmly with this lightweight, safety-adjust clamp. Holds 6.4 - 12.7 mm (1/4 - 1/2") diameter glass tubing/thermometers up to 178 mm (7") from support rod. Tightening wing nut applies tension to the nickel-plated phosphor bronze jaws. Overall length: 210 mm (8 1/4"). If attaching to rods, a 90° clamp holder is required. 1/pk