



Be Right™



Gas washing bottle, 1200 mL

Proizvod #:	2662200
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis—Hach Method 10065.

Specifikacije

Body Material:	Glass
Kapacitet:	1200 mL
Materijal:	Borosilicate Glass
Package size UOM:	1pk
Volumen/ veličina paketa:	1 each